


1000 Series: Instrument Kits

S1001 PI: Simpler Push in Instrument Kit

This kit contains the two hand instruments (Mallet & Seating Tool) needed to place Push in implants. It is used in conjunction with the S1003 ST: Simpler Threaded Surgical Kit which contains all the drills and other osteotomy hand instruments to place implants.

Contains:

- 1 - Mallet (S7012)
- 1 - Seating Tool (S7010)
- 1 - 3.25 mm Implant Retriever (S7008)
- 1 - 4.0 mm Implant Retriever (S7009)



S1002 IG: Simpler iGuide Instrument Kit

This kit contains all the drills and other hand instruments to determine the correct Location, Angulation, Diameter, Depth and Inter proximal distance between implants.

Contains:

- 2 - 1.5 mm Twist Drills (S6015)
- 2 - 2.0 mm Twist Drills (S6020)
- 1 - One Slot Centering Guide (S7041)
- 1 - Two Slot Centering Guide (S7022)
- 1 - Three Slot Centering Guide (S7043)
- 2 - 3.25 mm Alignment Jig (S7045)
- 2 - 4.0 mm Alignment Jig (S7046)
- 2 - 5.0 mm Alignment Jig (S7047)
- 2 - 6.0 mm Alignment Jig (S7048)
- 4 - 2.0 mm Depth Guide (S7020)
- 1 - Alignment Jig Holder (S7050)



S1003 ST: Simpler Threaded Surgical Instrument Kit

This kit contains all the osteotomy drills and other hand instruments to place implants.

Contains:

- 1 - 4.0 mm Tissue Punch for C/A (AS6007)
- 1 - Drill Extender (AS6005)
- 1 - 1.5 mm Twist Drills (S6015)
- 1 - 2.0 mm Twist Drills (S6020)
- 1 - 2.5 mm Twist Drills (S6025)
- 1 - 3.0 mm Tri Bladed Spade Drills (S6030-3)
- 1 - 3.25 mm Tri Bladed Spade Drills (S60325-3)
- 1 - 3.6 mm Tri Bladed Spade Drills (S6036-3)
- 1 - 4.0 mm Tri Bladed Spade Drills (S6040-3)
- 1 - 4.5 mm Tri Bladed Spade Drills (S6045-3)
- 1 - 5.0 mm Tri Bladed Spade Drills (S6050-3)
- 1 - 5.5 mm Tri Bladed Spade Drills (S6055-3) (Continued next page)



S1003 ST: Simpler Threaded Surgical Instrument Kit (Continued)

- 1 - 3.25 mm Tapered Internal Hex Driver for C/A (ST6010T)
- 1 - 4.0 mm Tapered Internal Hex Driver for C/A (ST6011T)
- 1 - Tapered 1.3 mm Alan Key (S7001T)
- 4 - Paralleling Pins (AS6015)
- 1 - Ratchet (S7014)
- 1 - 2.0 mm Depth Gauge (S7020 DG)
- 1 - 3.0 mm Depth Gauge (S7030 DG)
- 1 - 4.0 mm Depth Gauge (S7040 DG)
- 1 - 3.25 mm Internal Hex Driver for Ratchet (ST7030)
- 1 - 4.0 mm Internal Hex Driver for Ratchet (ST7031)

S1004 ND: Simpler Narrow Diameter Instrument Kit

- 1 - 1.5 mm Twist Drills (S6015)
- 1 - 2.0 mm Twist Drills (S6020)
- 4 - Paralleling Pins (AS6015)
- 1 - Ratchet (S7014)
- 1 - 2.0 mm Depth Gauge (S7020 DG)
- 1 - External Hex Driver for C/A (ST6012)
- 1 - External Hex Driver for Ratchet (ST7032U)

**S1005 R: Simpler Restorative Instrument Kit**

- 1 - External Hex Driver for Ratchet (ST7032)
- 1 - Hex Wrench (S7003)
- 1 - Tapered 1.3 mm Alan Key (S7001T)
- 1 - Torque Wrench (S7016)
- 1 - 1.3 mm Alan Key Driver for Ratchet (AS7033)

**S1006 D: Simpler Osteotomy Drill Kit**

- 1 - 4.0 mm Tissue Punch for C/A (AS6007)
- 1 - Drill Extender (AS6005)
- 1 - 1.5 mm Twist Drills (S6015)
- 1 - 2.0 mm Twist Drills (S6020)
- 1 - 2.5 mm Twist Drills (S6025)
- 1 - 3.0 mm Tri Bladed Spade Drills (S6030-3)
- 1 - 3.25 mm Tri Bladed Spade Drills (S60325-3)
- 1 - 3.6 mm Tri Bladed Spade Drills (S6036-3)
- 1 - 4.0 mm Tri Bladed Spade Drills (S6040-3)
- 1 - 4.5 mm Tri Bladed Spade Drills (S6045-3)
- 1 - 5.0 mm Tri Bladed Spade Drills (S6050-3)
- 1 - 5.5 mm Tri Bladed Spade Drills (S6055-3)

**1001 UND: Upgrade Components for Narrow Diameter Implants**

All the instruments and drills needed to upgrade the Surgical Threaded Kit (S1003 ST) to place Narrow Diameter Implants. The kit includes:

- 1 - C/A Driver for External Hex (ST6012)
- 1 - Ratchet Driver for External Hex (AT7032)